



REXLENE[®]

MONOFILAMENT POLYPROPYLENE SUTURES





















MONOFILAMENT POLYPROPYLENE SUTURES




















Polypropylene, a polymer of propylene
Non absorbable
High tensile strength
Resistance to repeated bending

These monofilament polypropylene sutures are blue colored for easy identification during surgery. polypropylene is a sterile monofilament non absorbable surgical suture material made from an isotactic crystalline stereoisomer of polypropylene. They also do not support infection as well as tissue rupture. Used for orthopaedic, plastic and micro surgeries, general closure and cardiovascular they also conform to the USP standards. Avoid damage to the surface of the suture.














Needle Type	Length & Color	Suture U.S.P (Metric)				
		10-0 (0.2)	9-0 (0.3)	8-0 (0.4)	7-0 (0.5)	6-0 (0.7)
3.7mm ⊙ 3/8	80MIC 15cm Blue	PPM1037				
4.5mm ⊙ 3/8	150MIC 15cm Blue		PPM945	PPM845		
4.7mm ⊙ 3/8	130MIC 15cm Blue		PPM947	PPM847		
5mm ⊙ 3/8	100MIC 15cm Blue	PPM1005	PPM905			
▽ 3/8	150MIC 30cm Blue	PPB1005D				
6mm ▽ 3/8	220MIC 15cm Blue		PPB906			
	220MIC 25cm Blue	PPB1006				
	220MIC 25cm Blue	PPB1006D				
	150MIC 30cm Blue	PPB1006T				
	150MIC 30cm Blue	PPB1006TD				
⊙ 3/8	90cm Blue		PPM906D	PPM806D	PPM706D	
6.5mm ⊙ 3/8	130MIC 15cm Blue			PPM865		
	45cm Blue			PPM865D		
7mm ⊙ 3/8	13cm Blue			PPM807S		
	45cm Blue			PPM807		
	45cm Blue				PPM707D	
	60cm Blue					PPM607D
8mm ⊙ 3/8	200MIC 15cm Blue			PPM808		
	220MIC 25cm Blue	PPM1008				
	220MIC 60cm Blue				PPM708PD	
	60cm Blue				PPM708D	PPM608D
	45cm Blue				PPM708SD	
	90cm Blue				PPM708LD	
9mm ▽ 3/8	220MIC 25cm Blue		PPE909			
	45cm Blue					PPE609
⊙ 3/8	45cm Blue				PPM709SD	
	60cm Blue				PPM709D	PPM609D
	60cm Blue					PPM609PD















POLYPROPYLENE

Needle Type	Length & Color	Suture U.S.P (Metric)					
		7-0 (0.5)	6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	
10mm  3/8 	45cm Blue	PPM710SD	PPM610SD	PPM510SD			
	50cm Blue	PPM710D					
	60cm Blue	PPM710LD	PPM610D	PPM510D			
	75cm Blue	PPM710BD					
1/2		45cm Blue	PPM710D				
11mm  3/8  3/8 	45cm Blue	PPE711	PPE611				
	45cm Blue		PPM611S				
	60cm Blue		PPM611				
	60cm Blue		PPM611D				
 1/2		45cm Blue		PPA611D			
12mm  3/8 	50cm Blue	PPM712	PPM612				
	75cm Blue		PPM612D	PPM512D			
13mm  3/8  3/8 	45cm Blue		PPE613	PPE513	PPE413		
	45cm Blue		PPM613S	PPM513S			
	45cm Blue			PPM513D			
	60cm Blue	PPM713D	PPM613D	PPM513D	PPM413D		
	75cm Blue		PPM613	PPM513			
	75cm Blue		PPM613LD	PPM513LD			
	1/2		60cm Blue	PPN713D	PPN613D		
	75cm Blue			PPN513D			
 3/8		45cm Blue			PPB413		
	50cm Blue			PPB513			
14mm  3/8 	50cm Blue			PPM514			
	75cm Blue				PPM414D		



















Needle Type	Length & Color	Suture U.S.P (Metric)				
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)
15mm  3/8  FINE 	45cm Blue	PPE615T	PPE515T	PPE415T		
	45cm Blue	PPE615	PPE515	PPE415	PPE315	
	75cm Blue				PPE315T	
 1/2 	90cm Blue		PPN515D	PPN415D	PPN315D	
16mm  3/8 	90cm Blue			PPM416D		
	45cm Blue		PPN516	PPN416		
17mm  1/2 	90cm Blue			PPT417	PPT317	
	90cm Blue					PPT217D
	75cm Blue					
 3/8 	75cm Blue			PPM417D		
	75cm Blue			PPM417		
 1/2 	75cm Blue		PPN517SD			
	90cm Blue	PPN617D	PPN517D	PPN417D		
18mm  3/8 	45cm Blue		PPE518	PPE418	PPE318	
	45cm Blue		PPE518T	PPE418T		
 1/2 	75cm Blue			PPN418		
	90cm Blue	PPN618D				
19mm  1/2 	90cm Blue			PPN419D		



















POLYPROPYLENE

Needle Type	Length & Color	Suture U.S.P (Metric)					
		6-0 (0.7)	5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	
20mm							
 1/2 	75cm Blue		PPN520D		PPN320D		
	90cm Blue			PPN420D			
	90cm Blue					PPN220	
21mm							
 3/8  <hr/> 1/2  <hr/>  3/8  <hr/>  1/2 	45cm Blue		PPE521	PPE421	PPE321		
	45cm Blue				PPF321		
	50cm Blue				PPB421		
	75cm Blue					PPN321	
22mm							
 3/8   <hr/> 1/2 	90cm Blue			PPM422D			
	100cm Blue					PPM222	
	75cm Blue		PPN522D	PPN422D	PPN322D		
	90cm Blue	PPN622D					










Needle Type	Length & Color	Suture U.S.P (Metric)				
		5-0 (1)	4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)
24mm						
 3/8	 90cm Blue		PPE424D			
 1/2	 90cm Blue			PPT324		
	 90cm Blue				PPT224D	
 1/2	 90cm Blue			PPN324D		
25mm						
 3/8	 45cm Blue		PPE425T	PPE325T		
	 45cm Blue	PPE525	PPE425	PPE325		
1/2	 70cm Blue			PPF325	PPF225	
 1/2	 75cm Blue		PPN425			
	 90cm Blue			PPN325D	PPN225D	

POLYPROPYLENE

Needle Type	Length & Color	Suture U.S.P (Metric)				
		4-0 (1.5)	3-0 (2)	2-0 (3)	0 (3.5)	1 (4)
25mm  1/2 	120cm Blue			PPN225		
26mm  1/2 	90cm Blue 90cm Blue		PPN426D PPN326D		PPN226	
28mm  1/2 	75cm Blue 90cm Blue			PPN328D PPN228D	PPN228	
29mm  3/8 	45cm Blue	PPE429				
30mm  1/2 	75cm Blue				PPN130	
31mm  1/2 	100cm Blue				PPN131	
33mm  1/2 	75cm Blue 100cm Blue		PPF333		PPF133	
35mm  1/2 	90cm Blue		PPN335D PPN235D			
36mm  1/2 	100cm Blue				PPT136	

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
37mm						
 3/8 	75cm Blue	PPE337				
1/2 	100cm Blue			PPF137		
40mm						
 1/2 	75cm Blue				PPN0140	
41mm						
 3/8 	100cm Blue	PPE241				
1/2 	100cm Blue			PPF141		
 1/2 	75cm Blue 100cm Blue				PPN0141	
					PPN141	
45mm						
 1/2 	75cm Blue				PPN0145	
46mm						
 1/2 	100cm Blue				PPT0146	
49mm						
 1/2 	100cm Blue				PPF0149	
55mm						
 	90cm Blue	PPS255				

POLYPROPYLENE

Needle Type	Length & Color	Suture U.S.P (Metric)				
		3-0 (2)	2-0 (3)	0 (3.5)	1 (4)	2 (5)
<div style="border: 1px solid black; padding: 2px; display: inline-block;">60mm</div>  	75cm Blue	PPS360	PPS260			
<div style="border: 1px solid black; padding: 2px; display: inline-block;">70mm</div>   	60cm Blue		PPS270D			
	75cm Blue		PPS270TD			
<div style="border: 1px solid black; padding: 2px; display: inline-block;">75mm</div>  	100cm Blue	PPS375	PPS275	PPS175		
<div style="border: 1px solid black; padding: 2px; display: inline-block;">90mm</div>  	90cm Blue		PPE290			

◆ LIGATURE AND NON-NEEDLE SUTURES

U.S.P (Metric)	POLYPROPYLENE (BLUE)		
	45cm × 12	60cm × 13	180cm × 1
4-0 (1.5)	PP-44512	PP-46013	PP-41801
3-0 (2)	PP-34512	PP-36013	PP-31801
2-0 (3)	PP-24512	PP-26013	PP-21801
1-0 (3.5)	PP-14512	PP-16013	PP-11801
1 (4)	PP-014512	PP-016013	PP-011801





SM Eng. Co., Ltd.
(주) 에스엠이엔지

356-3 SAMNAK-DONG SASANG-GU BUSAN KOREA TEL:81-51-305-8016 FAX:81-51-305-8021
E-mail : sales@sm-eng.net Website : www.sm-eng.net