



Rexpoint drilled end needles are made of stainless steel 300 series,

*** Advantage of 300 series REXPOINT NEEDLE**

1. Rexpoint 300 series needles are materialized a superior penetration performance which come from it's ideal geometric design with up-to-date coating technology.
2. Rexpoint 300 series needles are materialized easy and perfect attaching with unique conic works.
3. Rexpoint 300 series needles are materialized strong and unbreakable character with it's excellent ductility.



POINT TYPE

REVERSE CUTTING	CODE	CONVENTIONAL CUTTING	CODE	TAPER POINT	CODE
	F		E		N
TAPER CUTTING	CODE	TAPER CUTTING	CODE	BLUNT	CODE
	C		D		M
DIAMOND CUTTING	CODE	C-CUTTING	CODE	C-CUTTING	CODE
	J		G		H

SPATULA CUTTING	CODE	SPATULA CUTTING	CODE	PREMIUM SPATULAR	CODE
	A		B		I
PREMIUM CUTTING	CODE	PREMIUM CUTTING	CODE		
	K		L		

CIRCLE OF NEEDLE

1/4 CIRCLE	CODE	3/8 CIRCLE	CODE	7/16 CIRCLE	CODE
	Q		P		T
1/2 CIRCLE	CODE	5/8 CIRCLE	CODE	HOOK	CODE
	R		U		W
COMPOUND CIRCLE	CODE	STRAIGHT	CODE	SKI	CODE
	V		S		Y

DRILLED END NEEDLE ENCODING

F	080	36	40	R
①	②	③	④	⑤

- ① Needle shape : Reverse cutting
- ② Needle Dia (mm) : 0.80mm
- ③ Needle Length (mm) : 36mm
- ④ Needle Hole (mm) : 0.40 mm
- ⑤ Needle Circle : 1/2 circle

www.sm-eng.net

SM Eng

STAINLESS STEEL AISI 304, SUS 304
DRILLED END
SURGICAL NEEDLE

RexPoint

SM Eng Co., Ltd.
356-3 Samnak-Dong Sasang-Gu BUSAN KOREA
TEL: +82-51-305-8016 FAX: +82-51-305-8021
E-mail: sales@sm-eng.net, info@sm-eng.net

www.sm-eng.net

Needle attaching device.

Pneumatic type attaching machine has been developed and supply for our customer in order to attach our drilled needle to suture more easily and perfectly.
Each one dies and microscope with light are installed on attaching device.
Additional dies are optional.
In order to serve our customer, we SM Eng, also making and supplying other machines and tools/device related with suture manufacturing plant.