

SILICONE LOOPS

Product advantage

The Silicone Loops are used in surgical situations requiring retraction, identification or occlusion of arteries, veins, nerves, ureters or tendons. They are oval in cross-section, allowing for even pressure distribution, therefore minimizing damage to tissues.

The material used is platinum-cured, radiopaque, biocompatible silicone.

These items are soft and pliable, providing even pressure distribution, thus greatly minimizing trauma.

Technical Details

SIZE	MINI, MAXI		
RED	BLUE	YELLOW	WHITE
ARTERIES	VEINS	URETERS, TENDONS	NERVES

